

URGENT MEDICAL DEVICE RECALL

June 19, 2020

Product Name	Product Name
PowerLoc™ MAX Power-Injectable Infusion Set AllPoints™ Port Access System	SafeStep™ Huber Needle Set myPICK™ Port Access Kit
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	SafeStep™ Huber Needle Set Port Access Kit
PowerLoc™ MAX Power-Injectable Infusion Port Access Kit	SafeStep™ Huber Needle Set with GuardIva™ Antimicrobial Hemostatic Dressing
PowerLoc™ MAX Power-Injectable Infusion Set	SafeStep™ Huber Needle Set with GuardIva™ Antimicrobial Hemostatic Dressing Port Access Kit
PowerLoc™ MAX Power-Injectable Infusion Set with GuardIva™ Antimicrobial Hemostatic Dressing Port Access Kit	SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap
SafeStep™ Huber Needle Set	SafeStep™ Huber Needle Set Port Access Kit (PAK)
SafeStep™ Huber Needle Set AllPoints™ Port Access System	

Refer to Attachment A for All Catalog/Lot Numbers and Expiration Dates affected.

For the Attention of: Risk Manager, Materials Manager, Infection Control, Medical Director, Medical Device Safety Officer, Director of Nursing, Nursing Education, Pharmacy Directors

BD is initiating a voluntary recall for the PowerLoc™ MAX Power-Injectable Infusion Sets, SafeStep™ Huber Needle Sets, and the associated Port Access Kits containing the SafeStep™ Huber Needle Sets and the PowerLoc™ MAX Infusion Sets.

Description of the problem and health hazard(s):

BD has confirmed that the product catalog and lot numbers referenced in Attachment A may develop cracks or breaks in the tubing near the Luer™ or Y site of the device during use.

Devices that have developed cracks or breaks can cause leakage resulting in exposure to chemicals such as chemotherapeutics or biohazard material (e.g., blood). A crack or break could result in potential blood loss, catheter occlusion, air embolism, under or interrupted infusion, or site contamination which could lead to infection.

Our records indicate that you may have been shipped product referenced in this letter between February 2019 and May 2020.

PLEASE TAKE THE FOLLOWING ACTIONS:

1. Immediately review your inventory for the specific catalog and lot numbers referenced in Attachment A. Discontinue use of those devices and destroy all product subject to this recall following your institution's process for destruction.

2. Share this notification with all users of the product within your facility network to ensure they are also aware of this recall.
3. Complete the attached Customer Response Form and return to the BD contact noted on the form **whether or not you have any of the impacted material** so that BD may acknowledge your receipt of this notification per FDA requirements and subsequently process your credit accordingly for the discarded material.
 - a. Customers who purchased product directly from BD will receive credit.
 - b. Customers who did **not** purchase product directly from BD will need to contact their distributor to obtain credit.
4. Report any adverse health consequences experienced with the use of this product to BD at product.complaints@bd.com. Events may also be reported to the FDA's MedWatch Adverse Event Reporting program via:

Web: MedWatch website at www.fda.gov/medwatch

Phone: 1-800-FDA-1088 (1-800-332-1088)

Mail: MedWatch, HF-2, FDA, 5600 Fisher's Lane, Rockville, MD 20852-9787

Actions taken by BD:

BD has identified the root cause of this issue and initiated corrective actions to prevent recurrence of the issue.

Contact Information:

If you require further assistance, please contact:


BD Contact	US Contact Information
Customer/Technical Support	Phone: 800-290-1689 Monday – Friday between 8:00 am to 5:00 pm MDT SLC-BASD.FieldAction@bd.com

BD is committed to advancing the world of health. Our primary objectives are patient and user safety and providing you with quality products. We apologize for any inconvenience this issue may have caused you and thank you in advance for helping us to resolve this matter as quickly and effectively as possible.

Sincerely,



Klaus Hoerauf, MD
VP Global Medical Affairs
BD Medication Delivery Solutions



Gail Griffiths
Sr. Director Corporate Regulatory Compliance
BD US Region



Customer Recall Response Form
MDS-20-1979-FA
PowerLoc™ Max and SafeStep™

Please assist BD by promptly returning this form to BD Regulatory Compliance at:
 Email: BDRC8@bd.com or Fax No.: 1-312-949-0231

Facility: _____
Please use full, current facility name. Do not use initials.

Street Address: _____

City: _____ State: _____ Zip: _____

Recall Response Form Completed By:			
Account Number:		Telephone No.	
PO Number:		Fax No.	
Name:		Email:	
Title:			

Please check all that apply:

- I have read and understood the attached notice and have taken appropriate actions.
- We do not have any of the affected product(s) on hand.
- I have destroyed all affected product indicated below as available inventory at the time of receipt of this notification and request credit. Please indicate below how product was purchased:
 - I purchased product directly from BD
 - I purchased product through the following distributor and will work directly with my distributor to coordinate credit. (Distributor Name _____)

Product Name	Catalog No.	Lot No.	Quantity in Eaches

Please assist BD with assuring these communications are delivered to the appropriate person/function within your facility if that is not you.

Person/function responsible for the receipt and management of recalls at your facility:

Name: _____

Phone: _____

Email: _____

Attachment A - Affected Product List

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
PowerLoc™ MAX Power-Injectable Infusion Set	0131910	00801741047329	ASCXF007	10/24/2021	25/case
			ASCXF019	10/30/2021	25/case
			ASDRF074	4/22/2022	25/case
			ASDRF075	4/22/2022	25/case
			ASDRF079	4/24/2022	25/case
			ASDRF084	4/24/2022	25/case
			ASDVF021	8/10/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0131915	00801741047336	ASCXF001	10/21/2021	25/case
			ASDNF040	1/21/2022	25/case
			ASDNF043	1/22/2022	25/case
			ASDVF031	8/14/2022	25/case
			ASDWF118	9/27/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0131975	00801741047343	ASCXF006	10/23/2021	25/case
			ASCXF018	10/30/2021	25/case
			ASDRF053	4/16/2022	25/case
			ASDRF058	4/17/2022	25/case
			ASDRF060	4/17/2022	25/case
			ASDRF071	4/20/2022	25/case
			ASDVF029	8/13/2022	25/case
			ASDWF117	9/27/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0132010	00801741047350	ASCYF028	11/5/2021	25/case
			ASCYF029	11/6/2021	25/case
			ASCZF025	12/14/2021	25/case
			ASDPF029	2/24/2022	25/case
			ASDQF015	3/4/2022	25/case
			ASDQF039	3/16/2022	25/case
			ASDQF041	3/17/2022	25/case
			ASDRF065	4/18/2022	25/case
			ASDTF043	6/14/2022	25/case
			ASDTF054	6/16/2022	25/case
			ASDTF055	6/16/2022	25/case
			ASDTF061	6/17/2022	25/case
			ASDVF018	8/9/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0132015	00801741047367	ASCYF023	11/3/2021	25/case
			ASCYF024	11/3/2021	25/case
			ASDTF031	6/10/2022	25/case
			ASDTF032	6/10/2022	25/case
			ASDTF033	6/11/2022	25/case
			ASDTF037	6/12/2022	25/case
			ASDTF038	6/12/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0132075	00801741047374	ASCYF040	11/9/2021	25/case
			ASCYF076	11/13/2021	25/case
			ASCZF055	12/27/2021	25/case
			ASDNF001	1/1/2022	25/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			ASDNF002	1/1/2022	25/case
			ASDNF004	1/3/2022	25/case
			ASDNF018	1/8/2022	25/case
			ASDQF036	3/14/2022	25/case
			ASDTF075	6/19/2022	25/case
			ASDTF096	6/22/2022	25/case
			ASDTF104	6/25/2022	25/case
			ASDTF105	6/25/2022	25/case
			ASDTF115	6/28/2022	25/case
			ASDUF006	7/2/2022	25/case
			ASDUF019	7/5/2022	25/case
			ASDUF026	7/7/2022	25/case
			ASDUF033	7/9/2022	25/case
			ASDUF042	7/11/2022	25/case
			ASDUF059	7/14/2022	25/case
			ASDUF063	7/15/2022	25/case
			ASDUF078	7/18/2022	25/case
			ASDUF097	7/21/2022	25/case
			ASDUF102	7/23/2022	25/case
			ASDUF105	7/24/2022	25/case
			ASDUF111	7/27/2022	25/case
			ASDUF119	7/30/2022	25/case
			ASDVF005	8/2/2022	25/case
			ASDVF014	8/5/2022	25/case
			ASDWF070	9/13/2022	25/case
			ASDWF081	9/16/2022	25/case
			ASDWF086	9/19/2022	25/case
			ASDWF102	9/22/2022	25/case
			ASDWF114	9/26/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0132210	00801741047381	ASDTF023	6/9/2022	25/case
			ASDTF024	6/9/2022	25/case
			ASDTF025	6/9/2022	25/case
			ASDTF034	6/12/2022	25/case
			ASDTF035	6/12/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0132215	00801741047398	ASCYF020	11/3/2021	25/case
			ASDTF027	6/10/2022	25/case
			ASDTF028	6/10/2022	25/case
			ASDTF036	6/12/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0132275	00801741047404	ASCYF021	11/3/2021	25/case
			ASCYF022	11/3/2021	25/case
			ASDTF016	6/8/2022	25/case
			ASDTF017	6/8/2022	25/case
			ASDTF018	6/8/2022	25/case
			ASDTF019	6/8/2022	25/case
			ASDWF106	9/24/2022	25/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
PowerLoc™ MAX Power-Injectable Infusion Set	0141910	00801741047411	ASCXF009	10/25/2021	25/case
			ASCXF010	10/25/2021	25/case
			ASDRF088	4/25/2022	25/case
			ASDRF089	4/25/2022	25/case
			ASDRF092	4/26/2022	25/case
			ASDRF095	4/28/2022	25/case
			ASDRF096	4/28/2022	25/case
			ASDXF071	10/19/2022	25/case
			ASDXF073	10/20/2022	25/case
			ASENF043	1/10/2023	25/case
			ASENF047	1/11/2023	25/case
			ASENF048	1/11/2023	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0141915	00801741047428	ASCXF002	10/21/2021	25/case
			ASDNF042	1/21/2022	25/case
			ASDXF083	10/22/2022	25/case
			ASDXF084	10/22/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0141975	00801741047435	ASCXF005	10/22/2021	25/case
			ASCXF008	10/24/2021	25/case
			ASDSF028	5/13/2022	25/case
			ASDSF032	5/15/2022	25/case
			ASDSF044	5/17/2022	25/case
			ASDSF051	5/20/2022	25/case
			ASDSF055	5/22/2022	25/case
			ASDSF056	5/23/2022	25/case
			ASENF061	1/13/2023	25/case
			ASENF075	1/17/2023	25/case
			ASEQF074	3/20/2023	25/case
			ASEQF086	3/22/2023	25/case
ASEQF099	3/25/2023	25/case			
PowerLoc™ MAX Power-Injectable Infusion Set	0142010	00801741047442	ASDNF027	1/14/2022	25/case
			ASDNF031	1/16/2022	25/case
			ASDNF035	1/19/2022	25/case
			ASDNF038	1/20/2022	25/case
			ASDQF021	3/8/2022	25/case
			ASDQF025	3/10/2022	25/case
			ASDQF030	3/13/2022	25/case
			ASDVF088	8/29/2022	25/case
			ASDWF003	9/1/2022	25/case
			ASDWF009	9/2/2022	25/case
			ASDWF012	9/3/2022	25/case
			ASDWF019	9/6/2022	25/case
			ASDWF072	9/13/2022	25/case
			ASERF011	4/1/2023	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0142015	00801741047459	ASDQF002	2/28/2022	25/case
			ASDQF003	2/28/2022	25/case
			ASDQF004	2/28/2022	25/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			ASDQF006	3/3/2022	25/case
			ASDTF101	6/24/2022	25/case
			ASDVF084	8/28/2022	25/case
			ASDXF036	10/11/2022	25/case
			ASDXF037	10/11/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0142075	00801741047466	ASCZF056	12/27/2021	25/case
			ASCZF062	12/30/2021	25/case
			ASDNF003	1/2/2022	25/case
			ASDNF005	1/4/2022	25/case
			ASDNF029	1/15/2022	25/case
			ASDNF032	1/17/2022	25/case
			ASDNF037	1/20/2022	25/case
			ASDNF046	1/22/2022	25/case
			ASDNF052	1/25/2022	25/case
			ASDTF040	6/13/2022	25/case
			ASDTF044	6/15/2022	25/case
			ASDTF059	6/17/2022	25/case
			ASDTF076	6/20/2022	25/case
			ASDTF098	6/23/2022	25/case
			ASDTF106	6/25/2022	25/case
			ASDUF003	6/30/2022	25/case
			ASDUF009	7/3/2022	25/case
			ASDUF021	7/6/2022	25/case
			ASDUF028	7/8/2022	25/case
			ASDUF043	7/11/2022	25/case
			ASDUF060	7/14/2022	25/case
			ASDUF077	7/17/2022	25/case
			ASDUF095	7/20/2022	25/case
			ASDUF108	7/25/2022	25/case
			ASDUF112	7/27/2022	25/case
			ASDWF073	9/14/2022	25/case
			ASDWF082	9/17/2022	25/case
			ASDWF091	9/21/2022	25/case
			ASDWF092	9/21/2022	25/case
			ASDWF113	9/26/2022	25/case
			ASDWF120	9/27/2022	25/case
			ASDXF002	10/3/2022	25/case
			ASDXF015	10/8/2022	25/case
			ASDZF021	12/8/2022	25/case
			ASDZF029	12/11/2022	25/case
			ASDZF037	12/14/2022	25/case
			ASDZF047	12/17/2022	25/case
			ASENF007	1/4/2023	25/case
			ASENF069	1/15/2023	25/case
			ASENF079	1/19/2023	25/case
ASENF083	1/20/2023	25/case			
ASEPF007	2/16/2023	25/case			

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
PowerLoc™ MAX Power-Injectable Infusion Set	0142210	00801741047473	ASCYF018	11/3/2021	25/case
			ASCYF019	11/3/2021	25/case
			ASCZF026	12/14/2021	25/case
			ASDRF047	4/16/2022	25/case
			ASDRF048	4/16/2022	25/case
			ASDRF049	4/16/2022	25/case
			ASDRF050	4/16/2022	25/case
			ASDXF049	10/12/2022	25/case
			ASDXF050	10/12/2022	25/case
			ASENF026	1/8/2023	25/case
			ASENF027	1/8/2023	25/case
			ASENF028	1/8/2023	25/case
			ASERF006	1/4/2023	25/case
			ASERF007	1/4/2023	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0142215	00801741047480	ASCYF007	11/1/2021	25/case
			ASDRF063	4/18/2022	25/case
PowerLoc™ MAX Power-Injectable Infusion Set	0142275	00801741047497	ASCYF005	11/1/2021	25/case
			ASCYF006	11/1/2021	25/case
			ASDRF009	4/2/2022	25/case
			ASDRF015	4/5/2022	25/case
			ASDRF016	4/5/2022	25/case
			ASDRF028	4/8/2022	25/case
			ASDVF068	8/26/2022	25/case
			ASDZF060	12/20/2022	25/case
			ASDZF063	12/22/2022	25/case
			ASDZF066	12/22/2022	25/case
			ASDZF068	4/7/2023	25/case
			ASERF026	4/5/2023	25/case
			ASERF028	4/5/2023	25/case
			PowerLoc™ MAX Power-Injectable Infusion Set AllPoints™ Port Access System	AN132075	00801741118968
REDV3617	1/31/2021	10/case			
REDY1815	3/31/2021	10/case			
REDY2789	2/28/2021	10/case			
REDZ0385	3/31/2021	10/case			
REEN2192	6/30/2021	10/case			
REEN2782	6/30/2021	10/case			
PowerLoc™ MAX Power-Injectable Infusion Set AllPoints™ Port Access System	AN142010	00801741119231	REDW2682	2/28/2021	10/case
			REDX1471	3/31/2021	10/case
			REDX3857	3/31/2021	10/case
			REDY0765	4/30/2021	10/case
			REDY1593	4/30/2021	10/case
			REDY2943	4/30/2021	10/case
			REDZ0676	5/31/2021	10/case
			REDZ1903	5/31/2021	10/case
REDZ3077	5/31/2021	10/case			

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REEN1222	6/30/2021	10/case
			REEN3292	6/30/2021	10/case
			REEP0666	7/31/2021	10/case
			REEP4137	7/31/2021	10/case
			REER0612	9/30/2021	10/case
PowerLoc™ MAX Power-Injectable Infusion Set AllPoints™ Port Access System	AN142010G	00801741119323	REDW1021	11/30/2020	10/case
			REDX1360	11/30/2020	10/case
			REDX4728	2/28/2021	10/case
			REDY1816	4/30/2021	10/case
			REDY2935	4/30/2021	10/case
			REDY3765	4/30/2021	10/case
			REDZ1825	5/31/2021	10/case
			REDZ3146	5/31/2021	10/case
			REEN1935	2/28/2021	10/case
			REEP0773	6/30/2021	10/case
			REEP3544	7/31/2021	10/case
			REEP4066	6/30/2021	10/case
			REER0579	6/30/2021	10/case
			REER1103	9/30/2021	10/case
PowerLoc™ MAX Power-Injectable Infusion Set AllPoints™ Port Access System	AN142075	00801741119224	REDV2172	1/31/2021	10/case
			REDV3683	1/31/2021	10/case
			REDW0467	2/28/2021	10/case
			REDW0521	2/28/2021	10/case
			REDW1023	2/28/2021	10/case
			REDW1157	2/28/2021	10/case
			REDX1361	3/31/2021	10/case
			REDX1490	3/31/2021	10/case
			REDY0815	4/30/2021	10/case
			REDY1594	4/30/2021	10/case
			REDY2790	4/30/2021	10/case
			REDY3544	2/28/2021	10/case
			REDZ0110	5/31/2021	10/case
			REDZ0386	2/28/2021	10/case
			REDZ0677	5/31/2021	10/case
			REDZ1904	5/31/2021	10/case
			REDZ3078	5/31/2021	10/case
			REEN1072	6/30/2021	10/case
			REEN2186	6/30/2021	10/case
			REEN3075	6/30/2021	10/case
			REEN3310	6/30/2021	10/case
			REEP0663	7/31/2021	10/case
			REEP1401	7/31/2021	10/case
REEP4138	7/31/2021	10/case			
REEQ0611	8/31/2021	10/case			
REEQ1246	8/31/2021	10/case			

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REEQ2436	6/30/2021	10/case
			REEQ3589	6/30/2021	10/case
			REER0576	6/30/2021	10/case
			REER1046	9/30/2021	10/case
			REER2161	9/30/2021	10/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00001	00801741098345	REEN0212	6/30/2021	5/case
			REEP2685	7/31/2021	5/case
			REEP3660	7/31/2021	5/case
			REEQ3563	8/31/2021	5/case
			REER1968	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00002	00801741098352	REDZ0705	5/31/2021	5/case
			REEP3313	7/31/2021	5/case
			REEP3547	7/31/2021	5/case
			REER0577	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00028	00801741129728	REEN0214	6/30/2021	5/case
			REEQ3569	8/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00032	00801741137501	REDY3552	1/31/2021	5/case
			REEN1934	2/28/2021	5/case
			REEP3280	4/30/2021	5/case
			REEP3552	4/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00033	00801741137518	REER1115	7/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00036	00801741138911	REEN0218	6/30/2021	5/case
			REEP3659	7/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00051	00801741175718	REEP1436	6/30/2021	5/case
			REER0298	6/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00055	00801741182181	REER1335	5/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set myPICK™ Port Access Kit	CPA00056	00801741182907	REDY0992	4/30/2021	5/case
			REEN0278	6/30/2021	5/case
			REEP1437	4/30/2021	5/case
			REEQ0232	5/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2131910	00801741044151	REDW3798	2/28/2021	5/case
			REDY2866	4/30/2021	5/case
			REDZ1719	5/31/2021	5/case
			REEN3137	6/30/2021	5/case
			REEP0951	7/31/2021	5/case
			REEQ2684	8/31/2021	5/case
			REEQ4092	8/31/2021	5/case
			REER1229	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2131915	00801741044168	REEQ0664	8/31/2021	5/case
			REEQ1649	8/31/2021	5/case
			REER0459	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2131975	00801741044175	REDZ1818	5/31/2021	5/case
			REEP3049	7/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2132010	00801741044182	REDW1092	2/28/2021	5/case
			REDX3754	3/31/2021	5/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REDY1682	4/30/2021	5/case
			REDY2971	4/30/2021	5/case
			REDY3530	4/30/2021	5/case
			REDZ1070	5/31/2021	5/case
			REEN3271	6/30/2021	5/case
			REEN4480	6/30/2021	5/case
			REEP0992	7/31/2021	5/case
			REEP3957	7/31/2021	5/case
			REEQ1647	8/31/2021	5/case
			REEQ2724	8/31/2021	5/case
			REEQ4093	8/31/2021	5/case
			REER1266	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2132015	00801741044199	REDY3719	4/30/2021	5/case
			REEN1270	6/30/2021	5/case
			REEP0952	7/31/2021	5/case
			REEQ0684	8/31/2021	5/case
			REEQ2676	8/31/2021	5/case
			REEQ4054	8/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2132075	00801741044205	REEDT2879	11/30/2020	5/case
			REDU2256	12/31/2020	5/case
			REDU3371	12/31/2020	5/case
			REDU4487	12/31/2020	5/case
			REDV1088	1/31/2021	5/case
			REDV2972	1/31/2021	5/case
			REDV3470	1/31/2021	5/case
			REDW0402	2/28/2021	5/case
			REDW1087	2/28/2021	5/case
			REDX1518	3/31/2021	5/case
			REDY1804	4/30/2021	5/case
			REDY2867	4/30/2021	5/case
			REDY3531	4/30/2021	5/case
			REDZ1069	5/31/2021	5/case
			REEN1104	6/30/2021	5/case
			REEN4481	6/30/2021	5/case
			REEP0953	7/31/2021	5/case
			REEP3958	7/31/2021	5/case
			REEQ0665	8/31/2021	5/case
			REEQ1650	8/31/2021	5/case
REER0530	9/30/2021	5/case			
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2132210	00801741044212	REEN2139	6/30/2021	5/case
			REEP4053	7/31/2021	5/case
			REER0461	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2132275	00801741044236	REDY3720	4/30/2021	5/case
			REEQ4094	8/31/2021	5/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2141910	00801741044243	REDW3682	2/28/2021	5/case
			REDY0827	4/30/2021	5/case
			REDZ1071	5/31/2021	5/case
			REEP0994	7/31/2021	5/case
			REER0539	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2141915	00801741044250	REEN4551	6/30/2021	5/case
			REEP0954	7/31/2021	5/case
			REEQ1648	8/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2141975	00801741044267	REDW1197	2/28/2021	5/case
			REDY2972	4/30/2021	5/case
			REDZ3130	5/31/2021	5/case
			REEP0995	7/31/2021	5/case
			REEP2909	7/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2142010	00801741044274	REDW0546	2/28/2021	5/case
			REDX1519	3/31/2021	5/case
			REDX4930	3/31/2021	5/case
			REDY1683	4/30/2021	5/case
			REDY2868	4/30/2021	5/case
			REDY3532	4/30/2021	5/case
			REDZ1064	5/31/2021	5/case
			REDZ1819	5/31/2021	5/case
			REEN1046	6/30/2021	5/case
			REEN2146	6/30/2021	5/case
			REEP0955	7/31/2021	5/case
			REEP3052	7/31/2021	5/case
			REEP3959	7/31/2021	5/case
			REER0462	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2142015	00801741044281	REDZ1721	5/31/2021	5/case
			REEQ2678	8/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2142075	00801741044298	REDV1236	1/31/2021	5/case
			REDW0407	2/28/2021	5/case
			REDW0550	2/28/2021	5/case
			REDX4852	3/31/2021	5/case
			REDY0773	4/30/2021	5/case
			REDY0828	4/30/2021	5/case
			REDY1684	4/30/2021	5/case
			REDY2973	4/30/2021	5/case
			REDZ1063	5/31/2021	5/case
			REDZ1072	5/31/2021	5/case
			REDZ1820	5/31/2021	5/case
			REEN1271	6/30/2021	5/case
			REEN3273	6/30/2021	5/case
			REEP0956	7/31/2021	5/case
			REEP0996	7/31/2021	5/case
REEP1728	7/31/2021	5/case			

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REEP3053	7/31/2021	5/case
			REEP4054	7/31/2021	5/case
			REEQ1664	8/31/2021	5/case
			REEQ2720	8/31/2021	5/case
			REER0540	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2142210	00801741044304	REDW3683	2/28/2021	5/case
			REDY3721	4/30/2021	5/case
			REEN4552	6/30/2021	5/case
			REER1261	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2142215	00801741044311	REEN3155	6/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set Port Access Kit	2142275	00801741044328	REDW0403	2/28/2021	5/case
			REDY1685	4/30/2021	5/case
			REDZ1722	5/31/2021	5/case
			REEP0997	7/31/2021	5/case
			REEQ0666	8/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set with GuardIVa™ Antimicrobial Hemostatic Dressing Port Access Kit	2132010G	00801741095504	REDY0771	4/30/2021	5/case
			REDY2980	4/30/2021	5/case
			REEN1111	6/30/2021	5/case
			REEN4550	6/30/2021	5/case
			REER0529	9/30/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set with GuardIVa™ Antimicrobial Hemostatic Dressing Port Access Kit	2132015G	00801741095511	REEN2070	6/30/2021	5/case
			REEP0993	7/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set with GuardIVa™ Antimicrobial Hemostatic Dressing Port Access Kit	2141910G	00801741095542	REDX0120	11/30/2020	5/case
PowerLoc™ MAX Power-Injectable Infusion Set with GuardIVa™ Antimicrobial Hemostatic Dressing Port Access Kit	2142010G	00801741085680	REDW1093	11/30/2020	5/case
			REDW2530	11/30/2020	5/case
			REDX1388	3/31/2021	5/case
			REDY0772	4/30/2021	5/case
			REEN1112	6/30/2021	5/case
			REEN2779	2/28/2021	5/case
			REEP2910	7/31/2021	5/case
			REEQ4095	8/31/2021	5/case
PowerLoc™ MAX Power-Injectable Infusion Set with GuardIVa™ Antimicrobial Hemostatic Dressing Port Access Kit	2142015G	00801741085680	REER2358	9/30/2021	5/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
PowerLoc™ MAX Power-Injectable Infusion Set with GuardIva™ Antimicrobial Hemostatic Dressing Port Access Kit	2142210G	00801741085680	REDX0188	2/28/2021	5/case
SafeStep™ Huber Needle Set	LH-0029	00801741066160	ASCYF001	10/31/2021	25/case
			ASCYF017	11/3/2021	25/case
			ASDRF018	4/6/2022	25/case
			ASDRF031	4/9/2022	25/case
			ASDRF040	4/13/2022	25/case
			ASDSF065	5/27/2022	25/case
			ASDSF070	5/29/2022	25/case
			ASDTF004	6/1/2022	25/case
			ASDTF007	6/4/2022	25/case
			ASDVF048	8/19/2022	25/case
			ASDVF054	8/21/2022	25/case
			ASDXF058	10/14/2022	25/case
			ASDXF060	10/15/2022	25/case
			ASDXF066	10/17/2022	25/case
			ASDYF071	11/20/2022	25/case
			ASDYF077	11/23/2022	25/case
			ASDYF098	11/29/2022	25/case
			ASDZF010	12/3/2022	25/case
			ASDZF017	12/6/2022	25/case
			ASDZF024	12/9/2022	25/case
			ASDZF035	12/12/2022	25/case
			ASDZF036	12/13/2022	25/case
			ASDZF046	12/17/2022	25/case
			ASDZF071	12/24/2022	25/case
			ASDZF081	12/26/2022	25/case
			ASDZF093	12/30/2022	25/case
			ASENF002	1/2/2023	25/case
ASERF025	4/5/2023	25/case			
ASERF032	4/7/2023	25/case			
SafeStep™ Huber Needle Set	LH-0030	00801741066184	ASCYF015	11/2/2021	25/case
			ASCYF016	11/2/2021	25/case
			ASCZF027	12/14/2021	25/case
			ASDRF043	4/15/2022	25/case
			ASDRF044	4/15/2022	25/case
			ASDRF054	4/16/2022	25/case
			ASDRF055	4/17/2022	25/case
			ASDRF056	4/17/2022	25/case
			ASDVF058	8/23/2022	25/case
			ASDVF063	8/25/2022	25/case
			ASDXF020	10/9/2022	25/case
			ASDXF045	10/12/2022	25/case
			ASENF011	1/5/2023	25/case
			ASENF014	1/6/2023	25/case
ASERF004	4/1/2023	25/case			

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
SafeStep™ Huber Needle Set	LH-0031	00801741066207	ASCYF089	11/24/2021	25/case
			ASCYF090	11/25/2021	25/case
			ASCYF091	11/26/2021	25/case
			ASCYF092	11/26/2021	25/case
			ASCZF009	12/5/2021	25/case
			ASCZF010	12/5/2021	25/case
			ASCZF011	12/6/2021	25/case
			ASCZF022	12/12/2021	25/case
			ASCZF023	12/13/2021	25/case
			ASCZF039	12/18/2021	25/case
			ASCZF040	12/19/2021	25/case
			ASCZF045	12/22/2021	25/case
			ASCZF054	12/26/2021	25/case
			ASDPF016	2/19/2022	25/case
			ASDQF035	3/14/2022	25/case
			ASDQF054	3/22/2022	25/case
			ASDQF055	3/24/2022	25/case
			ASDQF060	3/27/2022	25/case
			ASDQF064	3/29/2022	25/case
			ASDRF006	4/1/2022	25/case
			ASDRF014	4/3/2022	25/case
			ASDRF019	4/7/2022	25/case
			ASDRF032	4/10/2022	25/case
			ASDRF039	4/13/2022	25/case
			ASDRF069	4/19/2022	25/case
			ASDRF073	4/22/2022	25/case
			ASDRF087	4/25/2022	25/case
			ASDRF098	4/28/2022	25/case
			ASDSF003	5/3/2022	25/case
			ASDSF010	5/6/2022	25/case
			ASDSF016	5/9/2022	25/case
			ASDSF019	5/10/2022	25/case
			ASDSF022	5/12/2022	25/case
			ASDSF033	5/16/2022	25/case
			ASDSF047	5/18/2022	25/case
			ASDSF050	5/19/2022	25/case
			ASDSF057	5/23/2022	25/case
			ASDSF061	5/26/2022	25/case
			ASDSF069	5/29/2022	25/case
			ASDTF005	6/1/2022	25/case
ASDTF008	6/4/2022	25/case			
ASDTF022	6/9/2022	25/case			
ASDVF051	8/20/2022	25/case			
ASDVF060	8/29/2022	25/case			
ASDVF071	8/26/2022	25/case			

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			ASDVF079	8/28/2022	25/case
			ASDVF091	8/29/2022	25/case
			ASDWF004	9/1/2022	25/case
			ASDWF013	9/4/2022	25/case
			ASDWF018	9/6/2022	25/case
			ASDWF034	9/10/2022	25/case
			ASDWF035	9/10/2022	25/case
			ASDWF036	9/10/2022	25/case
			ASDWF040	9/10/2022	25/case
			ASDWF041	9/10/2022	25/case
			ASDWF042	9/10/2022	25/case
			ASDWF046	9/11/2022	25/case
			ASDWF047	9/11/2022	25/case
			ASDWF048	9/11/2022	25/case
			ASDWF052	9/11/2022	25/case
			ASDWF053	9/11/2022	25/case
			ASDWF054	9/11/2022	25/case
			ASDWF058	9/11/2022	25/case
			ASDWF059	9/11/2022	25/case
			ASDWF060	9/11/2022	25/case
			ASDWF080	9/16/2022	25/case
			ASDWF090	9/20/2022	25/case
			ASDWF099	9/21/2022	25/case
			ASDWF121	9/27/2022	25/case
			ASDWF126	9/29/2022	25/case
			ASDXF006	10/5/2022	25/case
			ASDXF014	10/8/2022	25/case
			ASDXF053	10/13/2022	25/case
			ASDXF063	10/16/2022	25/case
			ASDXF070	10/19/2022	25/case
			ASDXF078	10/22/2022	25/case
			ASDXF092	10/24/2022	25/case
			ASDXF096	10/25/2022	25/case
			ASDXF101	10/31/2022	25/case
			ASDXF103	10/28/2022	25/case
			ASDXF106	10/29/2022	25/case
			ASDYF001	10/31/2022	25/case
			ASDYF004	11/1/2022	25/case
			ASDYF006	11/3/2022	25/case
			ASDYF009	11/4/2022	25/case
			ASDYF010	11/4/2022	25/case
			ASDYF015	11/6/2022	25/case
			ASDYF019	11/8/2022	25/case
			ASDYF023	11/9/2022	25/case
			ASDYF029	11/10/2022	25/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			ASDYF033	11/12/2022	25/case
			ASDYF037	11/13/2022	25/case
			ASDYF044	11/14/2022	25/case
			ASDYF050	11/16/2022	25/case
			ASDYF054	11/17/2022	25/case
			ASDYF060	11/18/2022	25/case
			ASDYF061	11/19/2022	25/case
			ASDYF074	11/21/2022	25/case
			ASDYF082	11/24/2022	25/case
			ASDYF093	11/27/2022	25/case
			ASDZF004	12/1/2022	25/case
			ASDZF013	12/4/2022	25/case
			ASDZF058	12/20/2022	25/case
			ASDZF069	12/23/2022	25/case
			ASEPF012	2/21/2023	25/case
			ASEPF015	2/23/2023	25/case
			ASEPF017	2/24/2023	25/case
			ASEPF021	2/26/2023	25/case
			ASEQF001	3/12/2023	25/case
			ASEQF022	3/6/2023	25/case
			ASEQF026	3/8/2023	25/case
			ASEQF031	3/9/2023	25/case
			ASEQF066	3/17/2023	25/case
			ASEQF083	3/21/2023	25/case
			ASEQF092	3/24/2023	25/case
			ASEQF110	3/27/2023	25/case
			ASEQF132	3/31/2023	25/case
			ASERF013	4/2/2023	25/case
			ASERF022	4/3/2023	25/case
			ASDNF034	1/19/2022	25/case
			ASDNF049	1/23/2022	25/case
			ASDNF055	1/27/2022	25/case
			ASDPF015	2/18/2022	25/case
			ASDPF017	2/19/2022	25/case
			ASDPF021	2/21/2022	25/case
			ASDPF024	2/23/2022	25/case
			ASDPF025	2/23/2022	25/case
			ASDPF026	2/23/2022	25/case
			ASDQF017	3/5/2022	25/case
			ASDRF066	4/18/2022	25/case
			ASDUF120	7/30/2022	25/case
			ASDVF004	8/1/2022	25/case
			ASDVF009	8/4/2022	25/case
			ASDVF013	8/5/2022	25/case
			ASDVF016	8/6/2022	25/case
SafeStep™ Huber Needle Set	LH-0032	00801741066221			

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			ASDVF025	8/11/2022	25/case
			ASDVF028	8/12/2022	25/case
			ASDVF032	8/14/2022	25/case
			ASDVF052	8/20/2022	25/case
			ASDXF048	10/12/2022	25/case
			ASDXF052	10/13/2022	25/case
			ASDXF061	10/15/2022	25/case
			ASDXF068	10/18/2022	25/case
			ASDXF077	10/21/2022	25/case
			ASDXF091	10/24/2022	25/case
			ASDXF100	10/27/2022	25/case
			ASDXF107	10/29/2022	25/case
			ASDYF083	11/25/2022	25/case
			ASDYF088	11/27/2022	25/case
			ASEPF008	2/18/2023	25/case
ASEPF011	2/20/2023	25/case			
ASEPF014	2/21/2023	25/case			
SafeStep™ Huber Needle Set	LH-0033	00801741066245	ASCXF012	10/26/2021	25/case
			ASCXF013	10/27/2021	25/case
			ASDSF006	5/3/2022	25/case
			ASDSF008	5/5/2022	25/case
			ASDSF012	5/7/2022	25/case
			ASDSF013	5/8/2022	25/case
			ASDSF021	5/11/2022	25/case
			ASDSF025	5/12/2022	25/case
			ASDWF110	9/25/2022	25/case
			ASDWF111	9/25/2022	25/case
			ASDWF116	9/26/2022	25/case
			ASDWF122	9/28/2022	25/case
			ASDXF003	10/4/2022	25/case
			ASENF058	1/13/2023	25/case
			ASEQF118	3/27/2023	25/case
SafeStep™ Huber Needle Set	LH-0034	00801741066269	ASCYF011	11/2/2021	25/case
			ASCYF012	11/2/2021	25/case
			ASDRF099	4/29/2022	25/case
			ASDRF100	4/29/2022	25/case
			ASDRF101	4/29/2022	25/case
			ASDSF001	5/2/2022	25/case
			ASDSF005	5/3/2022	25/case
			ASDXF010	10/7/2022	25/case
			ASDXF011	10/7/2022	25/case
			ASDXF016	10/8/2022	25/case
			ASDXF079	10/22/2022	25/case
			ASENF051	1/12/2023	25/case
SafeStep™ Huber Needle Set	LH-0035	00801741066283	ASCYF013	11/2/2021	25/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			ASCYF014	11/2/2021	25/case
			ASCZF029	12/14/2021	25/case
			ASDVF034	8/14/2022	25/case
			ASDVF035	8/15/2022	25/case
			ASDVF038	8/15/2022	25/case
			ASDVF042	8/16/2022	25/case
			ASDVF044	8/18/2022	25/case
			ASDXF054	10/13/2022	25/case
			ASDXF055	10/13/2022	25/case
			ASENF030	1/8/2023	25/case
			ASENF033	1/9/2023	25/case
			ASENF034	1/9/2023	25/case
			ASENF042	1/10/2023	25/case
SafeStep™ Huber Needle Set	LH-0036	00801741066306	ASCYF010	11/1/2021	25/case
			ASDRF059	4/17/2022	25/case
			ASDRF061	4/18/2022	25/case
			ASDRF062	4/18/2022	25/case
			ASDXF057	10/14/2022	25/case
			ASERF001	4/1/2023	25/case
			ASERF002	4/1/2023	25/case
			ASERF003	4/1/2023	25/case
SafeStep™ Huber Needle Set	LH-0037	00801741066320	ASDQF008	3/3/2022	25/case
			ASDQF009	3/3/2022	25/case
			ASDQF011	3/4/2022	25/case
			ASDQF012	4/3/2022	25/case
			ASDTF100	6/24/2022	25/case
			ASDVF075	8/27/2022	25/case
			ASDXF090	10/23/2022	25/case
			ASDXF093	10/24/2022	25/case
SafeStep™ Huber Needle Set	LH-0038	00801741066344	ASDPF043	2/27/2022	25/case
			ASDPF044	2/27/2022	25/case
			ASDPF045	2/27/2022	25/case
			ASDPF046	2/27/2022	25/case
			ASDQF001	2/28/2022	25/case
			ASDVF077	8/27/2022	25/case
			ASDVF078	8/28/2022	25/case
			ASDXF044	10/11/2022	25/case
			ASDXF047	10/12/2022	25/case
SafeStep™ Huber Needle Set	LH-0039	00801741066368	ASCXF003	10/21/2021	25/case
			ASDSF063	5/27/2022	25/case
			ASDXF080	10/22/2022	25/case
SafeStep™ Huber Needle Set	LH-0040	00801741066382	ASCXF004	10/21/2021	25/case
			ASDPF028	2/24/2022	25/case
			ASDSF064	5/27/2022	25/case
			ASDXF081	10/22/2022	25/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
SafeStep™ Huber Needle Set AllPoints™ Port Access System	AN152075	00801741119484	ASENF057	1/13/2023	25/case
			REDT0554	11/30/2020	10/case
			REDW1158	2/28/2021	10/case
			REDW2451	2/28/2021	10/case
			REDX1491	3/31/2021	10/case
			REDY1817	4/30/2021	10/case
			REDY2798	4/30/2021	10/case
			REEN3076	6/30/2021	10/case
			REEP3245	7/31/2021	10/case
			REEP3545	7/31/2021	10/case
SafeStep™ Huber Needle Set AllPoints™ Port Access System	AN152075G	00801741119583	REDR2659	9/30/2020	10/case
			REDS4631	10/31/2020	10/case
			REDT0709	11/30/2020	10/case
			REDT4233	11/30/2020	10/case
			REDU2308	12/31/2020	10/case
			REDV1300	11/30/2020	10/case
			REDV3619	1/31/2021	10/case
			REDX2864	11/30/2020	10/case
			REDY3766	2/28/2021	10/case
			REEP3542	7/31/2021	10/case
			REQQ0644	8/31/2021	10/case
			REER1047	6/30/2021	10/case
SafeStep™ Huber Needle Set AllPoints™ Port Access System	AN162010	00801741119682	REDX4729	12/31/2020	10/case
			REDZ0700	5/31/2021	10/case
			REDZ1905	5/31/2021	10/case
			REEN1900	6/30/2021	10/case
			REEN2797	6/30/2021	10/case
			REEP3240	7/31/2021	10/case
			REEP4079	7/31/2021	10/case
			REQQ0719	8/31/2021	10/case
			REQQ1251	8/31/2021	10/case
REER1105	9/30/2021	10/case			
SafeStep™ Huber Needle Set AllPoints™ Port Access System	AN162010G	00801741119446	REQQ1250	8/31/2021	10/case
SafeStep™ Huber Needle Set AllPoints™ Port Access System	AN162075	00801741119675	REDR2870	9/30/2020	10/case
			REDS0199	9/30/2020	10/case
			REDS3479	10/31/2020	10/case
			REDT1592	11/30/2020	10/case
			REDU0734	11/30/2020	10/case
			REDU1349	12/31/2020	10/case
			REDU2234	12/31/2020	10/case
			REDU3360	12/31/2020	10/case
			REDU4418	12/31/2020	10/case
			REDV1066	1/31/2021	10/case
			REDV3684	1/31/2021	10/case
			REDW0388	2/28/2021	10/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REDW1159	2/28/2021	10/case
			REDY1595	4/30/2021	10/case
			REDY2936	4/30/2021	10/case
			REDY3545	4/30/2021	10/case
			REDZ1827	5/31/2021	10/case
			REEN2784	2/28/2021	10/case
			REEP0777	6/30/2021	10/case
			REER0610	9/30/2021	10/case
			REER2163	9/30/2021	10/case
SafeStep™ Huber Needle Set AllPoints™ Port Access System	AN162075G	00801741119439	REDS4692	10/31/2020	10/case
			REDT0555	11/30/2020	10/case
			REDU1332	12/31/2020	10/case
			REDU2236	12/31/2020	10/case
			REDV1956	11/30/2020	10/case
			REDW3670	11/30/2020	10/case
			REDX1492	3/31/2021	10/case
			REDX3881	2/28/2021	10/case
			REDX4827	2/28/2021	10/case
			REDY1818	4/30/2021	10/case
			REDY2791	2/28/2021	10/case
			REDY3767	4/30/2021	10/case
			REDZ0388	2/28/2021	10/case
			REDZ1906	5/31/2021	10/case
			REDZ3147	5/31/2021	10/case
			REEN3311	6/30/2021	10/case
			REEQ1172	8/31/2021	10/case
			REER1106	9/30/2021	10/case
REER2125	9/30/2021	10/case			
SafeStep™ Huber Needle Set myPICK™ Port Access Kit	CPA00020A	00801741129827	REDZ0394	5/31/2021	5/case
			REDZ1664	5/31/2021	5/case
			REDZ3037	5/31/2021	5/case
			REEN1071	6/30/2021	5/case
			REER0578	9/30/2021	5/case
			REER1113	9/30/2021	5/case
SafeStep™ Huber Needle Set myPICK™ Port Access Kit	CPA00042A	00801741182969	REEN0220	6/30/2021	5/case
			REEN0271	6/30/2021	5/case
			REEN5242	6/30/2021	5/case
			REEP2688	7/31/2021	5/case
SafeStep™ Huber Needle Set myPICK™ Port Access Kit	CPA00049	00801741175695	REDZ1486	5/31/2021	5/case
			REDZ1678	5/31/2021	5/case
			REEP1382	6/30/2021	5/case
			REER1973	6/30/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0029	00801741085482	REDX2834	3/31/2021	5/case
			REDX4856	3/31/2021	5/case
			REDY1792	4/30/2021	5/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REDY2773	4/30/2021	5/case
			REEN2075	6/30/2021	5/case
			REEN4559	6/30/2021	5/case
			REEP1004	7/31/2021	5/case
			REEP3058	7/31/2021	5/case
			REEQ4021	8/31/2021	5/case
			REER0466	9/30/2021	5/case
			REER1235	9/30/2021	5/case
			REER2399	9/30/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0029YN	00801741085499	REDW3688	2/28/2021	5/case
			REDZ3026	5/31/2021	5/case
			REEN1109	6/30/2021	5/case
			REEP2916	7/31/2021	5/case
			REEQ2690	8/31/2021	5/case
			REEQ4063	8/31/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0030	00801741085505	REDS1399	10/31/2020	5/case
			REDS1875	10/31/2020	5/case
			REDS3371	10/31/2020	5/case
			REDS4425	10/31/2020	5/case
			REDT1836	11/30/2020	5/case
			REDU4512	12/31/2020	5/case
			REDX2861	3/31/2021	5/case
			REDY0778	4/30/2021	5/case
			REDY3542	1/23/2020	5/case
			REDZ3102	5/31/2021	5/case
			REEN1044	6/30/2021	5/case
			REEP3066	7/31/2021	5/case
			REEQ2648	8/31/2021	5/case
			SafeStep™ Huber Needle Set Port Access Kit	PA-0030YN	00801741085512
SafeStep™ Huber Needle Set Port Access Kit	PA-0031	00801741085529	REDR1897	9/30/2020	5/case
			REDR2695	9/30/2020	5/case
			REDR2796	9/30/2020	5/case
			REDS1288	10/31/2020	5/case
			REDS1400	10/31/2020	5/case
			REDS2408	10/31/2020	5/case
			REDS4621	10/31/2020	5/case
			REDT0697	11/30/2020	5/case
			REDT1645	11/30/2020	5/case
			REDT1837	11/30/2020	5/case
			REDT2924	11/30/2020	5/case
			REDT4207	11/30/2020	5/case
			REDU2262	12/31/2020	5/case
			REDU2347	12/31/2020	5/case
			REDU3328	12/31/2020	5/case
			REDV1094	1/31/2021	5/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REDV1239	1/31/2021	5/case
			REDV2772	1/31/2021	5/case
			REDV2983	1/31/2021	5/case
			REDW0413	2/28/2021	5/case
			REDX3752	3/31/2021	5/case
			REDY1723	4/30/2021	5/case
			REDY3731	4/30/2021	5/case
			REDZ0146	5/31/2021	5/case
			REDZ1859	5/31/2021	5/case
			REDZ3027	5/31/2021	5/case
			REEN2076	6/30/2021	5/case
			REEN2144	6/30/2021	5/case
			REEN3136	6/30/2021	5/case
			REEP0964	7/31/2021	5/case
			REER2366	9/30/2021	5/case
REER2367	9/30/2021	5/case			
REER2401	9/30/2021	5/case			
SafeStep™ Huber Needle Set Port Access Kit	PA-0031YN	00801741085536	REDN3937	6/30/2020	5/case
			REDQ0732	8/31/2020	5/case
			REDQ1532	8/31/2020	5/case
			REDQ3957	8/31/2020	5/case
			REDS1850	10/31/2020	5/case
			REDS2265	10/31/2020	5/case
			REDT0698	11/30/2020	5/case
			REDT1646	11/30/2020	5/case
			REDT2747	11/30/2020	5/case
			REDT4188	11/30/2020	5/case
			REDU0781	11/30/2020	5/case
			REDU2263	12/31/2020	5/case
			REDV2200	1/31/2021	5/case
			REDV3476	1/31/2021	5/case
			REDW0549	2/28/2021	5/case
			REDX1497	3/31/2021	5/case
			REDX3753	3/31/2021	5/case
			REDZ1056	5/31/2021	5/case
			REDZ1761	5/31/2021	5/case
			REEN1110	6/30/2021	5/case
REEN2794	6/30/2021	5/case			
REEN4560	6/30/2021	5/case			
REEP3956	7/31/2021	5/case			
SafeStep™ Huber Needle Set Port Access Kit	PA-0032	00801741085543	REDY3732	4/30/2021	5/case
			REDZ1860	5/31/2021	5/case
			REEN2774	6/30/2021	5/case
			REEN4491	6/30/2021	5/case
			REEP1005	7/31/2021	5/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REEQ1672	8/31/2021	5/case
			REER1236	9/30/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0032YN	00801741085550	REDQ1535	8/31/2020	5/case
			REDR2696	9/30/2020	5/case
			REDR2795	9/30/2020	5/case
			REDS1857	10/31/2020	5/case
			REDS4422	10/31/2020	5/case
			REDT1838	11/30/2020	5/case
			REDT2744	11/30/2020	5/case
			REDZ1762	5/31/2021	5/case
			REEN1045	6/30/2021	5/case
			REEN2786	6/30/2021	5/case
			REEN4561	6/30/2021	5/case
			REEP2919	7/31/2021	5/case
			REEP4051	7/31/2021	5/case
			REEQ4065	8/31/2021	5/case
			REER1273	9/30/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0033	00801741085567	REDX0146	3/31/2021	5/case
			REDY0779	4/30/2021	5/case
			REEN2792	6/30/2021	5/case
			REEP0965	7/31/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0033YN	00801741085574	REDZ3028	5/31/2021	5/case
			REEN4492	6/30/2021	5/case
			REEP0983	7/31/2021	5/case
			REEQ0701	8/31/2021	5/case
			REEQ1660	8/31/2021	5/case
			REEQ4024	8/31/2021	5/case
			REER0454	9/30/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0034	00801741085581	REDW1196	2/28/2021	5/case
			REDY0780	4/30/2021	5/case
			REDZ3104	5/31/2021	5/case
			REER2354	9/30/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0034YN	00801741085598	REDX1373	3/31/2021	5/case
			REEP0984	7/31/2021	5/case
			REEQ4066	8/31/2021	5/case
			REER2388	9/30/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0036	00801741085604	REEN1048	6/30/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0038	00801741085628	REEN2077	6/30/2021	5/case
			REEQ0671	8/31/2021	5/case
SafeStep™ Huber Needle Set Port Access Kit	PA-0038YN	00801741085635	REDX1374	3/31/2021	5/case
			REDZ1065	5/31/2021	5/case
			REEN4562	6/30/2021	5/case
			REEP0942	7/31/2021	5/case
SafeStep™ Huber Needle Set with GuardIva™ Antimicrobial Hemostatic Dressing	CPA00040	00801741143892	REEN0219	6/30/2021	5/case
			REEP1376	7/31/2021	5/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			REEP3657	6/30/2021	5/case
SafeStep™ Huber Needle Set with GuardIva™ Antimicrobial Hemostatic Dressing Port Access Kit	PA-0031G	00801741085666	REEN3159	6/30/2021	5/case
			REEP3067	7/31/2021	5/case
			REEQ4064	8/31/2021	5/case
SafeStep™ Huber Needle Set with GuardIva™ Antimicrobial Hemostatic Dressing Port Access Kit	PA-0031YNG	00801741085673	REDS3541	10/31/2020	5/case
			REDV3438	11/30/2020	5/case
			REEN2781	6/30/2021	5/case
SafeStep™ Huber Needle Set with GuardIva™ Antimicrobial Hemostatic Dressing Port Access Kit	PA-0032G	00801741087639	REEN2145	6/30/2021	5/case
			REEP3068	7/31/2021	5/case
SafeStep™ Huber Needle Set with GuardIva™ Antimicrobial Hemostatic Dressing Port Access Kit	PA-0033G	00801741096433	REDY0835	4/30/2021	5/case
			REDZ3103	5/31/2021	5/case
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0029YN	00801741066177	ASDQF068	3/30/2022	25/case
			ASDRF001	3/31/2022	25/case
			ASDRF002	3/31/2022	25/case
			ASDRF003	3/31/2022	25/case
			ASDRF005	4/1/2022	25/case
			ASEPF018	2/25/2023	25/case
			ASEPF019	2/25/2023	25/case
			ASEPF020	2/26/2023	25/case
			ASEPF026	3/11/2023	25/case
			ASEQF008	3/4/2023	25/case
			ASEQF025	3/7/2023	25/case
			ASEQF030	3/9/2023	25/case
			ASEQF039	3/13/2023	25/case
			ASEQF041	4/13/2023	25/case
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0030YN	00801741066191	ASCYF038	11/8/2021	25/case
			ASCYF097	11/28/2021	25/case
			ASDQF065	3/30/2022	25/case
			ASDQF066	3/30/2022	25/case
			ASDQF067	3/30/2022	25/case
			ASDWF061	9/10/2022	25/case
			ASDWF062	9/10/2022	25/case
			ASDWF063	9/10/2022	25/case
			ASDWF064	9/10/2022	25/case
			ASDWF065	9/10/2022	25/case
			ASEQF045	3/14/2023	25/case
			ASEQF046	3/14/2023	25/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0031YN	00801741066214	ASCYF039	11/9/2021	25/case
			ASCYF065	11/12/2021	25/case
			ASCYF093	11/27/2021	25/case
			ASCZF007	12/3/2021	25/case
			ASCZF008	12/4/2021	25/case
			ASCZF018	12/9/2021	25/case
			ASCZF019	12/9/2021	25/case
			ASCZF041	12/20/2021	25/case
			ASCZF042	12/21/2021	25/case
			ASDNF026	1/13/2022	25/case
			ASDQF027	3/11/2022	25/case
			ASDQF029	3/12/2022	25/case
			ASDQF037	3/16/2022	25/case
			ASDQF043	3/19/2022	25/case
			ASDQF051	3/21/2022	25/case
			ASDQF061	3/27/2022	25/case
			ASDRF037	4/12/2022	25/case
			ASDYF011	11/4/2022	25/case
			ASDYF016	11/7/2022	25/case
			ASDYF024	11/9/2022	25/case
			ASDYF034	11/12/2022	25/case
			ASDYF045	11/15/2022	25/case
			ASDYF055	11/17/2022	25/case
			ASDYF064	11/19/2022	25/case
			ASDYF066	11/20/2022	25/case
			ASDYF075	11/22/2022	25/case
			ASDYF078	12/4/2022	25/case
			ASDYF087	11/26/2022	25/case
			ASDZF073	12/24/2022	25/case
			ASDZF077	12/25/2022	25/case
			ASDZF087	12/30/2022	25/case
			ASDZF094	12/12/2022	25/case
			ASENF004	1/2/2023	25/case
			ASENF009	1/5/2023	25/case
			ASENF020	1/7/2023	25/case
			ASENF040	1/10/2023	25/case
			ASENF056	1/19/2023	25/case
			ASENF071	1/16/2023	25/case
			ASENF077	1/19/2023	25/case
			ASENF078	1/19/2023	25/case
ASENF086	1/22/2023	25/case			
ASENF087	1/25/2023	25/case			
ASENF088	1/27/2023	25/case			
ASENF089	1/29/2023	25/case			
ASEPF001	2/4/2023	25/case			

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
			ASEPF002	2/6/2023	25/case
			ASEPF003	2/9/2023	25/case
			ASEPF004	2/11/2023	25/case
			ASEPF005	2/13/2023	25/case
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0032YN	00801741066238	ASCYF025	11/4/2021	25/case
			ASCYF026	11/4/2021	25/case
			ASCYF027	11/5/2021	25/case
			ASCZF028	12/14/2021	25/case
			ASDPF030	2/25/2022	25/case
			ASDQF016	3/5/2022	25/case
			ASDQF019	3/7/2022	25/case
			ASDQF024	3/10/2022	25/case
			ASDRF036	4/12/2022	25/case
			ASDZF014	12/4/2022	25/case
			ASDZF019	12/6/2022	25/case
			ASDZF023	12/9/2022	25/case
			ASDZF034	12/12/2022	25/case
			ASDZF040	12/14/2022	25/case
			ASDZF044	12/16/2022	25/case
			ASDZF052	12/18/2022	25/case
			ASDZF059	12/20/2022	25/case
			ASDZF067	12/22/2022	25/case
			ASEQF056	3/16/2023	25/case
			ASEQF069	3/18/2023	25/case
ASEQF075	3/20/2023	25/case			
ASEQF076	3/20/2023	25/case			
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0033YN	00801741066252	ASCXF016	10/28/2021	25/case
			ASCXF017	10/29/2021	25/case
			ASDQF007	3/3/2022	25/case
			ASDQF010	3/3/2022	25/case
			ASDQF013	3/4/2022	25/case
			ASERF021	4/23/2023	25/case
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0034YN	00801741066276	ASCXF014	10/28/2021	25/case
			ASCXF015	10/28/2021	25/case
			ASDPF038	2/26/2022	25/case
			ASDPF039	2/27/2022	25/case
			ASEQF108	3/26/2023	25/case
			ASEQF109	3/26/2023	25/case
			ASEQF111	3/27/2023	25/case
			ASEQF117	3/27/2023	25/case
			ASEQF121	3/28/2023	25/case
			ASEQF126	3/30/2023	25/case
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0035YN	00801741066290	ASCYF096	11/28/2021	25/case
			ASDWF083	9/17/2022	25/case
			ASDWF084	9/17/2022	25/case

Product Name	Catalog Number	UDI (GTIN, DI+PI)	Lot Number	Expiration Date (MM/DD/YYYY)	Package Size
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0036YN	00801741066313	ASEQF044	3/14/2023	25/case
			ASDRF007	4/2/2022	25/case
			ASDRF008	4/2/2022	25/case
			ASEQF042	3/14/2023	25/case
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0037YN	00801741066337	ASCYF008	11/1/2021	25/case
			ASCYF009	11/1/2021	25/case
			ASCYF095	11/28/2021	25/case
			ASCZF030	12/14/2021	25/case
			ASDPF031	2/25/2022	25/case
			ASDPF032	2/25/2022	25/case
			ASCYF003	11/1/2021	25/case
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0038YN	00801741066351	ASCYF004	11/1/2021	25/case
			ASDPF036	2/26/2022	25/case
			ASDPF037	2/26/2022	25/case
			ASEQF085	3/22/2023	25/case
			ASEQF093	3/24/2023	25/case
			ASEQF094	3/24/2023	25/case
			ASEQF095	3/24/2023	25/case
			ASERF041	4/29/2023	25/case
			SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0039YN	00801741066313
ASCYF094	11/28/2021	25/case			
ASDPF040	2/27/2022	25/case			
ASDPF041	2/27/2022	25/case			
SafeStep™ Huber Needle Set with Y-Injection Site and Medegen™ Needleless Injection Cap	LH-0040YN	00801741066399	ASCXF011	10/26/2021	25/case
			ASDNF044	1/22/2022	25/case
			ASDNF045	1/22/2022	25/case
			ASEQF100	3/25/2023	25/case
			ASEQF101	3/25/2023	25/case
			ASEQF102	3/25/2023	25/case
			ASEQF103	3/25/2023	25/case
			ASEQF104	3/25/2023	25/case